

# COPPER ION SENSOR

## Authors

Dr hab. Joanna Makowska  
Dr hab. Krzysztof Żamojć  
Dr hab. Dariusz Wyrzykowski  
Dominik Kamrowski  
Prof. Wiesław Wiczak  
Prof. Lech Chmurzyński

Department of Chemistry  
University of Gdańsk

## Commercialization opportunities



- ➔ Granting of the licenece
- ➔ Transfer of ownership
- ➔ Joint further research or commercialization

## IP Status



The invention was submitted for patenting according to Polish (P.430245) procedures.

## Implementation progress



TRL 7

Copper ions ( $\text{Cu}^{2+}$  ions) - the third most abundant transition metal in the human body - play an important role in many fundamental physiological processes.

Both the insufficiency of  $\text{Cu}^{2+}$  ions, which is, among others, a symptom of anaemia, as well as their excess, are harmful and can lead to many diseases, such as Alzheimer's, Parkinson's, Huntington's diseases or those related to the digestive system, liver or blood vessels. Due to the wide use of copper, which is a heavy metal, e.g. in electrical engineering, electronics, construction and motorization or fertilizers, copper ions are also an important cause of environmental pollution.

The scientists from University of Gdańsk developed an innovative sensor (a new fluorescent peptide) for determination of copper ions in solutions.

### Advantages:

- technology does not entail the use of expensive or specialised equipment,
- very high sensitivity: possibility of detecting even single molecules,
- very high selectivity: a peptide, tested for interactions with approximately 20 metal cations, gives a response exclusively under the influence of copper ions.

### Main applications:

- analysing the quality of chemical reagents used in laboratories,
- measuring the concentration of copper ions: in waters (including mineral and deionised waters), in food drinks, in cosmetic and pharmaceutical products as well as in biological material for diagnostic use.

## Technology Transfer Office



tto@ug.edu.pl



+48 58 523 33 74

+48 58 523 33 75



ul. Jana Bażyńskiego 1a  
80-309 Gdańsk